

Art Robbins Instruments is pleased to announce the release of the new 3-well Intelli-Plate for protein crystallization. This plate features three 1.0µL reaction wells and a 140µL screen reservoir. Each identical well features a round bottom for easy crystal harvesting.

Contact us or our distributors for samples or for more information.

Hampton Research Catalog Number: HR3-183 - Per Pack: 10
HR3-185 - Per Case: 40

Hampton Research
34 Journey
Aliso Viejo, CA 92656-3317 U.S.A.
Tel: (949) 425-1321 • Fax: (949) 425-1611
Website: www.hamptonresearch.com

New 3-well Intelli-Plate

